

Application Note

MagXtract® 3200, a fully automated nucleic acid extraction and PCR setup system using MACHEREY-NAGEL's NucleoMag® Blood 200µl kit.



CE IVD FDA



Magnetic Bead

Fast and efficient automated nucleic acid purifications from various bio-specimens.



PCR Setup

Automatically and precisely dispense eluates, PCR reagents into PCR plate.



Open Platform

Comprehensive protocols creation and adaption for purification and PCR setup.



Whirl Mixing

Patented whirl-and-stir mixing approach to maximize the efficiency of purifications.




UV / HEPA

Safeguard against the contamination with UV light and HEPA.



Cool Station

Cooling control for eluates storage and reagent kits to maintain samples stability.



NucleoMag® Blood 200 µl

Automated extraction of genomic DNA from human whole blood on the MagXtract 3200

Introduction

MACHEREY-NAGEL developed the magnetic bead based NucleoMag® Blood 200µl kit allowing the automated purification of genomic DNA from whole blood. Chroma ATE provides high quality automated instrument, MagXtract 3200, for nucleic acid extraction. Through using NucleoMag® Blood 200µl combining with MagXtract 3200, laboratory scientists will conduct clinical examinations in a safe, accurate, and efficient procedure.

MagXtract 3200 Specification

Description	Automated nucleic acid extraction and liquid handling system
Capacity	Max. 32, 16-sample per run
Weight	70 kg (± 5kg)
Dimensions	721mm x 530mm x 567mm (W x D x H)

Excellent Performance

We evaluated the performance of purification of genomic DNA from whole blood using the NucleoMag® Blood 200µl kit and MagXtract 3200. High yield and purity of extracted

genomic DNA from whole blood are shown below.







Figure. High yield and purity of extracted genomic DNA from whole blood. A human whole blood sample (ZenBio, =W36982100056100) was used as a sample for the NucleoMag® Blood 200µl kit for 16 repeats and followed the extraction procedure on the MagXtract 3200 instrument. Genomic DNA were measured by Nanodrop ONE (ThermoFisher). A. Concentration, B. Yield, C, Ratio 260/280, D, Ratio 260/230


Open System and Automation Solution

MagXtract 3200 allows the authorized user to edit protocols including nucleic acid extraction and PCR preparation. High flexibility application can fit most magnetic based extraction kits and PCR reagent kits.

Open System – Based on MACHEREY-NAGEL extraction kit IFU, optimize the detailed factor to make excellent performance.


Automation Solution – Replace the manual steps including not only extraction but also PCR preparation.

The layout of extraction buffer in 96 deep well plate as below picture.



Please make sure the column 1 and 7 of deep well plate are empty. Because run 1 spin tip will be discarded in these column after lysis and then MagXtract 3200 will automatically add binding buffer and magnetic beads into column 2 and 8. Each specimen need two spin tips. Please make sure your spin tips are sufficient.

We have customized Elute and Storage Deck for 20 or 30ml bottle to aliquot binding buffer. The layout of reagents as below GUI picture.



Fully automatic steps:


1. Adding specimen and proteinase K to lysis buffer (column 2 and 8)
2. Specimen lysis
3. Adding binding buffer and magnetic beads
4. Collecting beads
5. Washing
6. Vaporing (drying)

7. Eluting
8. Transferring eluate to storage vial
9. Dispensing qPCR premix to each vial
10. Adding template to each vial

It saves about 90% hands-on time!

Semi Automation Solution – For customer who purchased standard elute and storage deck, we provide an alternative method to conduct semi-auto extraction.

The layout of extraction buffer in 96 deep well plate as below picture.



Please make sure the column 1 and 7 of deep well plate are empty. Because run 1 spin tip will be discarded in these column after lysis and then pause automatically to let user add binding buffer and magnetic beads.

Each specimen need two spin tips. Please make sure your spin tips are sufficient.

Semi automation steps:

1. Adding specimen to lysis buffer (column 2 and 8)
2. Specimen lysis
3. Automatically pause
4. **Manually** added binding buffer and magnetic beads and remember to remove first run spin tip from column 1 / 7.
5. Collecting beads
6. Washing
7. Vaporizing (drying)
8. Eluting
9. Transferring eluate to storage vial
10. Dispensing qPCR premix to each vial
11. Adding template to each vial

MagXtract 3200 GUI

Three operating modes: full run, preparation and extraction.

MagXtract 3200 software provides the protocol-based control to streamline the workflows. The stepwise GUI and touchscreen control guide the user through the complication of the assay setup, from sample loading to the consumable placement.

For more information please contact

MACHEREY-NAGEL Bioanalysis technical

support: bio-tech@mn-net.com

Chroma ATE Technical support:

tech-support.LSI@chroma.com.tw

Order Information

Product	REF
NucleoMag® Blood	744501.1 (96 preps)
200µl	744501.4 (384 preps)
MagXtract 3200	9-49200002
(CE-IVD)	
Spin tip in box	A94-000135
Filtered tip (50µl)	A94-000130
Filtered tip (1000µl)	A94-000131
96 deep well plate	Please contact Chroma ATE
20ml bottle	Please contact Chroma ATE
30ml bottle	Please contact Chroma ATE
Customized accessories	
Sample drawer	13 x 75 mm / 16 x 100 mm
Elute and storage deck	We can modify the number of hole to fit your application