



## **Application Note**

MagXtract® 3200, a fully automated nucleic acid extraction and PCR setup system using MACHEREY-NAGEL's NucleoMag® DNA Microbiome kit.







### **Magnetic Bead**

Fast and efficient automated nucleic acid purifications from various bio-specimens.



#### **Whirl Mixing**

Patented whirl-and-stir mixing approach to maximize the efficiency of purifications.



## **PCR Setup**

Automatically and precisely dispense eluates, PCR reagents into PCR plate.



#### **UV / HEPA**

Safeguard against the contamination with UV light and HEPA.




## **Open Platform**

Comprehensive protocols creation and adaption for purification and PCR setup.



#### **Cool Station**

Cooling control for eluates storage and reagent kits to maintain samples stability.









# NucleoMag® DNA Microbiome

**Dimensions**  $721 \text{mm} \times 530 \text{mm} \times 567 \text{mm}$  (W x D x H)

#### Automated extraction of bacteria DNA

# from stool sample on the **MagXtract**3200

#### Introduction

The commensal microbiome is important for human and animal health. Several studies indicate that specific changes in the microbiome involve in the development of various diseases. Molecular biological examination can provide useful data to understand microbiome. MACHEREY-NAGEL developed the magnetic bead based NucleoMag® DNA Microbiome kit allowing the automated purification of nucleic acids from stool or soil specimen. Chroma ATE provides high quality automated instrument, MagXtract 3200, for nucleic acid extraction. Through using NucleoMag® DNA Microbiome kit combining with MagXtract 3200, laboratory scientists will conduct examinations in a safe, accurate, and efficient procedure.

#### **MagXtract 3200 Specification**

Description	Automated nucleic acid	
	extraction and liquid	
	handling system	
Capacity	Max. 32, 16-sample per run	
Weight	70 kg (± 5kg)	

#### **Excellent Performance**

We evaluated the performance of purification of genomic DNA from 50mg stool samples using the NucleoMag® DNA Microbiome kit and MagXtract 3200. High yield and purity of extracted genomic DNA from stools are shown below table. The extracted DNA for downstream qPCR application is also conducted as below figure.


Sample	Conc.	Yield	Purity	Purity
(50mg)	(ng/µl)	(µg)	(260/280)	(260/230)
Human	19.3	1.6	1.92	1.63
stool				
Dog	41.6	3.3	1.97	2.09
stool				
Cat	16.1	1.3	1.98	1.59
stool				

Table. High yield and purity of extracted genomic DNA from stools. Different source of stool were added into NucleoMag® DNA Microbiome kit and followed the extraction procedure on the MagXtract 3200 instrument. Genomic DNA quality were evaluated using Nanodrop ONE (ThermoFisher).






**Figure. Downstream qPCR application.** All the gDNA template in qPCR reaction are normalized to 0.5, 5 and 50ng/reaction. And these extracted samples are detected *Bifidobacterium* positive using qPCR. A. Human stool, B. Dog stool, C. Cat stool.








#### **Open System and Automation Solution**

MagXtract 3200 allows the authorized user to edit protocols including nucleic acid extraction and PCR preparation. High flexibility application can fit most magnetic based extraction kits and PCR reagent kits.

**Open System** – Based on MACHEREY-NAGEL extraction kit IFU, optimize the detailed factor to make excellent performance.

**Automation Solution** – Replace the manual steps including not only extraction but also PCR preparation.

Note. According to NucleoMag® DNA Microbiome IFU, the lysates should be processed manually in advance. The layout of extraction buffer in 96 deep well plate for automation, please refer to below picture.









1. Binding

2. Collecting beads

3. Washing x 4

4. Vaporing (drying)

5. Eluting

6. Transferring eluate to storage vial

7. Dispensing qPCR premix to each vial

8. Adding template to each vial

It saves about 90% hands-on time!

#### MagXtract 3200 GUI

Three operating modes: full run, preparation and extraction.

MagXtract 3200 software provides the protocolbased control to streamline the workflows. The stepwise GUI and touchscreen control guide the user through the complication of the assay setup, from sample loading to the consumable placement.

For more information please contact

MACHEREY-NAGEL Bioanalysis technical

support: bio-tech@mn-net.com

**Chroma ATE** Technical support: tech-support.LSI@chroma.com.tw

#### **Order Information**

Product	REF	
NucleoMag <sup>®</sup> DNA	744330.1 (96 preps)	
Microbiome	744330.4 (384 preps)	
MagXtract 3200	9-49200002	
(CE-IVD)		
Spin tip in box	A94-000135	
Filtered tip (50µl)	A94-000130	
Filtered tip (1000µl)	A94-000131	
96 deep well plate	Please contact Chroma ATE	
Customized accessories		
Sample drawer	13 x 75 mm / 16 x 100 mm	
Elute and storage	We can modify the number	
deck	of hole to fit your application	

