



Application Note

MagXtract[®] 3200, a fully automated nucleic acid extraction and PCR setup system using MACHEREY-NAGEL's NucleoMag[®] Virus kit.





Magnetic Bead

Fast and efficient automated nucleic acid purifications from various bio-specimens.



Whirl Mixing

Patented whirl-and-stir mixing approach to maximize the efficiency of purifications.



PCR Setup

Automatically and precisely dispense eluates, PCR reagents into PCR plate.



UV / HEPA

Safeguard against the contamination with UV light and HEPA.



Open Platform


CE IVD FD/A

Comprehensive protocols creation and adaption for purification and PCR setup.



Cool Station

Cooling control for eluates storage and reagent kits to maintain samples stability.









NucleoMag® Virus

Automated extraction of SARS-CoV-2
RNA from swab matrix on the

MagXtract 3200

Introduction

Owing to COVID-19 pandemic, the demands of nucleic acid purification elevate than before and the needs of laboratory automation increase dramatically. MACHEREY-NAGEL developed the magnetic bead based NucleoMag® Virus kit allowing the automated purification of nucleic acids from specimen. Chroma ATE provides high quality automated instrument, MagXtract 3200, for nucleic acid extraction. Through using NucleoMag® Virus kit combining with MagXtract 3200, laboratory scientists will conduct clinical examinations in a safe, accurate, and efficient procedure.

MagXtract 3200 Specification

Description	Automated nucleic acid extraction and liquid handling system
Capacity	Max. 32, 16-sample per run
Weight	70 kg (± 5kg)
Dimensions	721mm x 530mm x 567mm (W x D x H)

Excellent Performance

We evaluated the performance of purification of SARS-CoV-2 viral RNA from swab matrix sample using the NucleoMag[®] Virus kit and MagXtract 3200. Highly sensitive target signal was shown as figure.




Figure. Highly sensitive detection of SARS-CoV-2 RNA in swab matrix. A serial dilution of inactivated SARS-CoV-2 virus (Amplirun® total SARS CoV-2 Control (Swab), Vircell microbiologists) was added into NucleoMag® Virus kit and followed the extraction procedure on the MagXtract 3200 instrument. Viral RNA was detected using qRT-PCR kit (Direct SARS-CoV-2 Realtime PCR kit, Vircell microbiologists) on the thermalcycler, LightCycler 96, Roche (Amplirun 1:20 = 1500 copies/ml, 15 copies/qPCR vial;







Amplirun 1:80 = 375 copies/ml, 3.75 copies/qPCR vial).


Open System and Automation Solution

MagXtract 3200 allows the authorized user to edit protocols including nucleic acid extraction and PCR preparation. High flexibility application can fit most magnetic based extraction kits and PCR reagent kits.

Open System – Based on MACHEREY-NAGEL extraction kit IFU, optimize the detailed factor to make excellent performance.

Automation Solution – Replace the manual steps including not only extraction but also PCR preparation.

The layout of extraction buffer in 96 deep well plate as below picture.




Please make sure the column 1 and 7 of deep well

plate are empty. Because run 1 spin tip will be discarded in these column after lysis and then MagXtract 3200 will automatically add binding buffer and magnetic beads into column 2 and 8.

Each specimen need two spin tips. Please make sure your spin tips are sufficient.

We have customized Elute and Storage Deck for 20 or 30ml bottle to aliquot binding buffer. The layout of reagents as below GUI picture.



Fully automatic steps:

- 1. Adding specimen and <u>proteinase K and carrier RNA</u> to lysis buffer (<u>note. proteinase K and carrier RNA must premix in advance to place on elute and storage deck)</u>
- 2. Specimen lysis
- 3. Adding binding buffer and magnetic beads
- 4. Collecting beads
- 5. Washing
- 6. Vaporing (drying)
- 7. Eluting
- 8. Transferring eluate to storage vial
- 9. Dispensing qPCR premix to each vial
- 10. Adding template to each vial

It saves about 90% hands-on time!








Semi Automation Solution – For customer who purchased standard elute and storage deck, we provide an alternative method to conduct semiauto extraction.

The layout of extraction buffer in 96 deep well plate as below picture.



Please make sure the column 1 and 7 of deep well plate are empty. Because run 1 spin tip will be

discarded in these column after lysis and then pause automatically to let user add binding buffer

and magnetic beads.

Each specimen need two spin tips. Please make sure your spin tips are sufficient.

Semi automation steps:

- Adding specimen to lysis buffer (column 2 and
 adding specimen to lysis buffer (column 2 and
- 2. Specimen lysis
- 3. Automatically pause

- 4. *Manually* added binding buffer and magnetic beads and <u>remember to remove first run spin tip</u> from column 1 / 7.
- 5. Collecting beads
- 6. Washing
- 7. Vaporing (drying)
- 8. Eluting
- 9. Transferring eluate to storage vial
- 10. Dispensing qPCR premix to each vial
- 11. Adding template to each vial

MagXtract 3200 GUI

Three operating modes: full run, preparation and extraction.

MagXtract 3200 software provides the protocolbased control to streamline the workflows. The stepwise GUI and touchscreen control guide the user through the complication of the assay setup, from sample loading to the consumable placement.

For more information please contact

MACHEREY-NAGEL Bioanalysis technical

support: bio-tech@mn-net.com

Chroma ATE Technical support: tech-support.LSI@chroma.com.tw







Order Information

Product	REF
NucleoMag [®] Virus	744800.1 (96 preps)
	744800.4 (384 preps)
MagXtract 3200	9-49200002
(CE-IVD)	
Spin tip in box	A94-000135
Filtered tip (50µl)	A94-000130
Filtered tip (1000µl)	A94-000131
96 deep well plate	Please contact Chroma ATE
20ml bottle	Please contact Chroma ATE
30ml bottle	Please contact Chroma ATE
Customized accessories	
Sample drawer	13 x 75 mm / 16 x 100 mm
Elute and storage	We can modify the number
deck	of hole to fit your application

